



Zirconia translucent ST/HT

Indications

- Partial / fully anatomical crown and bridge up to 16 units with maximum 2 pontics
- Crown and bridge up to 16 units for veneering
- Implant abutments
- Telescopic work and attachment
- Implant-supported superstructure

Technical details	Parameters (bridge)
<ul style="list-style-type: none"> • Flexural strength: 1200 – 1470 Mpa • Light Transmittance: 30-40 % • Shrinkage Factor: 19-22 % 	<ul style="list-style-type: none"> • Margin thickness: 0,25 mm • Wall thickness: 0,60 mm • Connector section: 9 - 15 mm²

Colour

A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4

! The part in the blank that is responsible for stability is in the dentine-shaded area. Since 80% of the bridge connectors should be positioned there, depending on the situation, there may be less incisal parts and thus an overall darker result in the restoration.

Our tip: Consider to select one shade lighter when working with a pronounced Spee curve, connectors at the level of the occlusal plane and restorations with a low overall height (up to 8 mm)

MILLING

ARTM ST/HT Block (Super/High Translucent)	Preis	KATANA Zirconia Block	Preis
Coping	€ 32,90	Coping	€ 50,20
Full anatomie	€ 32,90	Full anatomie	€ 50,20
Inlay/Onlay	€ 32,90	Inlay/Onlay	€ 50,20
Veneer	€ 32,90	Veneer	€ 50,20
Bridge	€ 32,90	Bridge	€ 50,20
Implant		Implant	
Coping	€ 32,90	Coping	€ 58,20
Full anatomie	€ 32,90	Full anatomie	€ 62,20

Scan / Design service

On the basis of the model you send in, which we scan for you, or on the basis of the transmitted scan data (intraoral scan or model scan), our experienced CAD specialists create your desired design and implement it in the appropriate material.

Coping	4,10 €
Vollanatomie	6,20 €
Inlay/Onlay	3.50 €
Veneer	3.50 €
Bridge	6,00 €
Standby	-2,00 €
OneDay	+6,20 €
Express Service - Goodwill insurance	+4,10 €

Standby	Requirement
Standby means that you benefit from an unbeatable price and leave us an additional working day of production time.	for all materials (except partial framework / tertiary structure) <ul style="list-style-type: none"> only by delivery of design file, intraoral scan or model scan

One Day	Requirement
Selecting this option your order will be dispatched the same day, if the following requirements are met:	Milling zirconia (ARTM translucent + CADtools SHT Color) <ul style="list-style-type: none"> Single coping / bridge up to 3 units Deadline of files Mo.-Fr. by 9:00 a.m.

Express Service Milling

For urgent orders which arrive after the deadline and shall be produced the same day. However, prior consultation by phone is required.

Goodwill insurance	Requirement
<p>Insure your work against impression errors and inconveniences with this all-round carefree package: in the event of damage, we will produce the affected structure again.</p> <p>The goodwill insurance applies in the following cases:</p> <ul style="list-style-type: none"> • - mistake in the impression from your dentist • - design errors • - wrong color choice • - accidental damage in your laboratory, etc. 	<ul style="list-style-type: none"> • only by delivery of file • the design lies within the indication • the work is newly commissioned • the newly commissioned work is in the same material <ul style="list-style-type: none"> • the goodwill notification takes place within one year, for splints within one month • you send us the affected work back <p>The credit note will be issued as soon as we have the damaged work in our center. Errors by shipping service providers and special services which have been ordered and carried out without warranty are excluded.</p> <p>Refunds via goodwill insurance are up to max. 3.5 % of total net sales possible (sales from delivery of files for the respective financial year).</p>

Production times / Data delivery

Milling Zirconia	Deadline			Shipment ²
	Mo. - Th.	Fr.	Sa. ¹	
Crown / bridge / etc. up to 5 units	4 p.m.	2 p.m.	11 a.m.	next working day
Crown / bridge / etc. from 6 units and in peak hours	4 p.m.	2 p.m.	11 a.m.	second working day
Zirconia custom shaded	4 p.m.	2 p.m.	-	second working day

¹ Only correct STL files will be produced, as there is no telephone customer service.

² Standard delivery time; except peak / rush hours. In this case a later delivery may occur.

SHIPMENT

Shipment outside Austria ¹	Delivery	Cost
DHL Express - End of Day ca. 18:00 h	Mo. - Fr.	depending on destination country
DHL Express - 13:00 h ²	Mo. - Fr.	
DHL Saturday delivery ²	Sa.	

¹ Please consider that delivery may occur in 24 - 72 hours, according on the country and region your laboratory is located and custom procedures, if needed.

² Option available only in some areas of the national territory. Contact our customer service team to check if the option is available in your municipality.

Payment options / Further information

Payment options

- Term of payment :
5 (five) working days after the date of shipment, exclusively via SEPA direct debit mandate.
- Grant us a direct debit authorization and you will receive 2 % cash discount on the invoice amount. Debiting will be effected every 10th of the following month.
- From € 4,000.00 net sales per month, we grant you a 5 % volume discount.

ARTAMED Milling Center

ADRESSE: Barbara-Prammer-Allee 15, 1220, Wien, Österreich
 Tel: +43 660 986 5820 / 01 9972747
 E-Mail: Office@artamed.at