



MEDIT

i900

ALL IN

Comfort
Control
Care

Everything you need to take
your practice to the next level.

Visit [getmedit.com/ALL-IN](https://www.getmedit.com/ALL-IN) or scan
the QR code to learn more.



i9000

ALL IN

Comfort • Control • Care

Discover what makes the Medit i900 THE BEST IN ITS CLASS.

Experience superior scanning technology with the new Medit i900. Featuring:

- faster, smoother scanning
- easier handling
- improved scan accuracy

Our next-level tech gives you more precision and control while the sleek design provides ultimate comfort and ease of use. Plus, with our unparalleled support and no hidden fees, **the i900 helps you focus on delivering outstanding patient care.**

INTUITIVE CONTROL

Real-time haptic and light feedback, 360° TouchBand™ navigation, and customizable command actions via our interactive touch pad. The i900 is our smartest, most accessible scanner yet.

ADVANCED ERGONOMICS

The i900 is smaller, lighter, and symmetrically-designed so you can scan from any position with absolute comfort.

SUPERIOR VISUAL OUTPUT

A 10-bit camera, wider FOV, and our 3rd Generation Optical Engine give you outstanding scan accuracy and resolution, along with errorless metal and edentulous scanning.

HASSLE FREE

Jointless unibody design eliminates bio-hazard concerns, and a light, flexible cable provides more freedom of movement. Your i900 also comes with a pre-configured laptop, reducing service issues and ensuring optimal scanner performance.



Grow your practice efficiently with i900's unparalleled technology.

The Medit i900 and our Medit Link suite of applications create an unmatched combination of performance and potential. No hidden fees or subscriptions - just one system with everything needed for streamlined workflows and patient care. It includes an optimized laptop, support, and a 3-year warranty, delivering security and solutions all-in-one package.



Get started with the i900 today.

Visit getmedit.com/ALL-IN
or scan the QR code to learn more.

MEDIT